

SURGI-TEC



THE ART OF MOVING BONE



WWW.SURGI-TEC.COM

SURGI-TEC

THE ART OF MOVING BONE



Surgi-Tec is an ISO 13485 certified medical device manufacturer based in Belgium. We specialise in CMF surgery implants.

The development of the TPD (transpalatal distractor) and the OBA (ortho bone anchor) have caused a remarkable shift in the field of CMF surgery.

These products are developed in close collaboration with Prof. Dr. M.Y. Mommaerts. This made the company renowned in Belgium and its neighbouring countries.

Quality is our top priority: we only work with certified suppliers for all our materials.

Our high-precision machines and the expertise of our production team assures that these quality standards are guaranteed in our assembly process.

The same quality requirements are continued in the relationship with our partners and customer service. Personal contact is one of the core values of Surgi-Tec. Our innovation and know-how have led to a very extensive OEM and subcontracting business.

Your idea is our common project!

SURGI-TEC NV BELGIUM
POORTAKKERSTRAAT 43-9051 SINT-DENIJS-WESTREM
PHONE: +32 (0)9 230 05 69
WWW.SURGI-TEC.COM - INFO@SURGI-TEC.COM

TABLE OF CONTENTS

TRANS PALATAL DISTRACTORS

TPD NEO	7
TPD ALL-IN-ONE	9
TPD CLASSIC	11
INSTRUMENTS AND SCREWS TPD	12

TRANS MANDIBULAR DISTRACTORS

TMD	16
INSTRUMENTS AND SCREWS TMD	17

ORTHODONTIC BONE ANCHORS

MOMMAERTS	19
MAXILLA	21
MANDIBULA	22
INSTRUMENTS AND SCREWS OBA	23

OSTEOSYNTHESIS PLATES

OSP	26
INSTRUMENTS AND SCREWS OSP	27

OSTEOSYNTHESIS SCREWS

OS	30
----	----


INSTRUMENTS

DRILLS	33
MOMMAERTS	34
SURGI-TEC	36

DISTRIBUTION

MEDARTIS	38
INFINITAS	39
META	40

TRANS PALATAL DISTRACTOR NEO



EXTRA SHARP
SPIKES

COUNTER NUT

CLOSED BODY DESIGN

COLOR MARKING
| STEP IS 0.33 mm DISTRACTION



INDICATIONS FOR
PLACEMENT

3 DIMENSIONAL MOVEMENT

ROUNDED EDGES

SLIM AND CLOSED DESIGN INCREASES PATIENT COMFORT

SURGI-TEC 

THE ART OF MOVING BONE



SURGI-TEC TRANS PALATAL DISTRACTOR


TPD NEO

TITANIUM GRADE 2 & 5

	REF.	DESCRIPTION	INITIAL WIDTH	DISTRACTION	MAX. WIDTH
	03-701A	TPD NEO MODULE 1	17 mm	7 mm	24 mm
	03-702A	TPD NEO MODULE 2	20 mm	12 mm	32 mm
	03-725A	TPD NEO MODULE 2.5	24 mm	19 mm	43 mm
	03-703A	TPD NEO MODULE 3	27 mm	25 mm	52 mm

TRANS PALATAL DISTRACTOR

ALL-IN-ONE



EXTRA SHARP
SPIKES

LOCKING SCREW



INDICATIONS
FOR PLACEMENT

COLOR MARKING
1 STEP IS 0.33 mm DISTRACTION

ROUNDED EDGES

3 DIMENSIONAL MOVEMENT



TELESCOPIC DESIGN FOR
MAXIMUM DISTRACTION

SURGI-TEC

THE ART OF MOVING BONE




SURGI-TEC TRANS PALATAL DISTRACTOR

TPD ALL-IN-ONE

TITANIUM GRADE 2 & 5


	REF.	DESCRIPTION	INITIAL WIDTH	DISTRACTION	MAX. WIDTH
	03-901A	TPD AIO MODULE 1	15 mm	7 mm	22 mm
	03-902A	TPD AIO MODULE 2	18 mm	12 mm	30 mm
	03-925A	TPD AIO MODULE 2.5	21 mm	18 mm	39 mm
	03-903A	TPD AIO MODULE 3	24 mm	23 mm	47 mm
	03-904A	TPD AIO MODULE 4	29 mm	29 mm	58 mm

TRANS PALATAL DISTRACTOR CLASSIC



EXTRA SHARP
SPIKES

LOCKING SCREW



SMOOTH SURFACE

COLOR MARKING
1 STEP IS 0.33 mm DISTRACTION

ROUNDED EDGES

HOLES FOR LIGATURE

3 PIECE DESIGN ALLOWS FOR INTERCHANGEABLE MODULES

SURGI-TEC 
THE ART OF MOVING BONE



SURGI-TEC TRANS PALATAL DISTRACTOR

TPD CLASSIC

TITANIUM GRADE 2 & 5

	REF.	DESCRIPTION	INITIAL WIDTH	DISTRACTION	MAX. WIDTH
	03-800A	ABUTMENT PLATES		2.5 mm holes	
	03-801A	TPD CLASSIC MODULE 1	15 mm	7 mm	22 mm
	03-802A	TPD CLASSIC MODULE 2	18 mm	12 mm	30 mm
	03-825A	TPD CLASSIC MODULE 2.5	21 mm	18 mm	39 mm
	03-803A	TPD CLASSIC MODULE 3	24 mm	23 mm	47 mm
	03-804A	TPD CLASSIC MODULE 4	29 mm	29 mm	58 mm



SURGI-TEC TRANS PALATAL DISTRACTOR

INSTRUMENTS

TPD NEO

	REF.	DESCRIPTION
	03-751A	TPD NEO - HINGED KEY
	03-750S	TPD NEO - PATIENT KEY
	99-901A	HANDPIECE FOR ALL INSERTS
	99-909S	INSERT - INTERNAL PENTAGON
	99-910S	INSERT - EXTERNAL PENTAGON



DRILLS CAN BE FOUND ON PAGE 33 OF THE CATALOGUE



SURGI-TEC TRANS PALATAL DISTRACTOR

INSTRUMENTS

TPD ALL-IN-ONE AND TPD CLASSIC

	REF.	DESCRIPTION
	03-951A	TPD ALL-IN-ONE / CLASSIC HINGED KEY
	03-950S	TPD ALL-IN-ONE / CLASSIC PATIENT KEY
	99-901A	HANDPIECE FOR ALL INSERTS
	99-101A	SCREWDRIVER- INTERNAL PENTAGON
	99-909S	INSERT - INTERNAL PENTAGON
	99-910S	INSERT - EXTERNAL PENTAGON



DRILLS CAN BE FOUND ON PAGE 33 OF THE CATALOGUE



SURGI-TEC TRANS PALATAL DISTRACTOR

INSTRUMENTS

DUMMIES (SCALE 1:1)

	REF.	DESCRIPTION
	03-790S	TPD NEO DUMMIES
	03-990S	TPD ALL-IN-ONE / CLASSIC DUMMIES

SCREWS

TO USE WITH TPD

* ALL TPD'S ARE COMPATIBLE WITH 2.5 mm RESCUE SCREWS

	REF.	DIMENSIONS	HEAD TYPE
TPD NEO (2 screws needed)			
 	70-707S	Ø 2.0 mm - LENGTH 7 mm	PENTAGON
TPD ALL-IN-ONE (2 screws & 1 locking screw needed)			
 	70-707S	Ø 2.0 mm - LENGTH 7 mm	PENTAGON
 	99-100S	LOCKING SCREW	INTERNAL PENTAGON
TPD CLASSIC (4 screws & 1 locking screw needed)			
 	70-707S	Ø 2.0 mm - LENGTH 7 mm	PENTAGON
 	99-100S	LOCKING SCREW	INTERNAL PENTAGON



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides.
If you prefer to receive our screws in 5-packed format, please add .5 to the reference.
e.g.: 70-205S in 5-pack becomes 70-205S.5

TRANS MANDIBULAR DISTRACTOR

DISTRACTION POSSIBLE ON 2 SIDES

0.5 mm PITCH

2 CARDANS ALLOWS
FOR SYMMETRIC
AND / OR
ASYMMETRIC
DISTRACTION

37.5°



FLAT

ANGLED HOLES

47.5°



BULKY

SURGI-TEC

THE ART OF MOVING BONE



SURGI-TEC TRANS MANDIBULAR DISTRACTOR

TMD

TITANIUM GRADE 2 & 5

	REF.	DESCRIPTION	DISTRACTION
	06-131A	TMD PEDIATRIC SHORT MODEL	12 mm
	06-132A	TMD PEDIATRIC LONG MODEL	20 mm
	06-141A	TMD BULKY SHORT MODEL	12 mm
	06-142A	TMD BULKY LONG MODEL	20 mm
	06-241A	TMD FLAT SHORT MODEL	12 mm
	06-242A	TMD FLAT LONG MODEL	20 mm



SURGI-TEC TRANS MANDIBULAR DISTRACTOR

INSTRUMENTS

TMD

	REF.	DESCRIPTION
	99-901A	HANDPIECE FOR ALL INSERTS
	99-909S	INSERT - INTERNAL PENTAGON
	99-910S	INSERT - EXTERNAL PENTAGON
	06-950A	TMD PATIENT KEY

SCREWS

RECOMMENDED TO USE WITH TMD

	REF.	DIMENSIONS	HEAD TYPE
TMD (6 screws needed)			
	70-809S	Ø 2.3 mm - LENGTH 9 mm	PENTAGON
	70-811S	Ø 2.3 mm - LENGTH 11 mm	PENTAGON
	70-813S	Ø 2.3 mm - LENGTH 13 mm	PENTAGON




All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5



DRILLS CAN BE FOUND ON PAGE 33 OF THE CATALOGUE

ORTHODONTIC BONE ANCHOR MOMMAERTS



SLA: SAND
BLASTED &
ACID ETCHED

TRIANGLE DESIGN FOR ZYGOMA

Ø 2.4 MM

ROUNDED EDGES

0.6 MM PLATE THICKNESS

EXTRA STRONG DESIGN
BUT REMAINS EASY BENDABLE
INTRA OPERATIVE

2 ROUNDED HOOKS

WIRE HOLE

BRACKET
SLOT: 0.56 x 0.76 mm / 22 x 30
TORQUE: -1°, ANGULATION: 0°

TUBES
SLOT: 2 x 0.56 x 0.70 mm / 22 x 28
TORQUE: -10°, DISTAL OFFSET: 0°

2 TYPES

SURGI-TEC

THE ART OF MOVING BONE



OBA MOMMAERTS

TITANIUM GRADE 2

	REF. NO.	DESCRIPTION	NECK
HOOKS ONLY			
	02-003S	OBA SLA - SHORT NECK - DOUBLE HOOK	4,5 mm
	02-006S	OBA SLA - LONG NECK - DOUBLE HOOK	6,5 mm
	02-009S	OBA SLA - ZYGOMA - DOUBLE HOOK	12,5 mm

02-003S 02-006S 02-009S

HOOKS AND BRACKET

	02-002A	OBA SLA - SHORT NECK - BRACKET	4,5 mm
	02-005A	OBA SLA - LONG NECK - BRACKET	6,5 mm
	02-008A	OBA SLA - ZYGOMA - BRACKET	12,5 mm


02-002A 02-005A 02-008A

HOOKS AND TUBE

	02-001A	OBA SLA - SHORT NECK - TUBE	4,5 mm
	02-004A	OBA SLA - LONG NECK - TUBE	6,5 mm
	02-007A	OBA SLA - ZYGOMA - TUBE	12,5 mm

02-001A 02-004A 02-007A

ORTHODONTIC BONE ANCHOR





SURGI-TEC ORTHODONTIC BONE ANCHOR

OBA MAXILLA

TITANIUM GRADE 2

REF.	DESCRIPTION	NECK
 02-005S	02-005S LONG NECK - DOUBLE HOOK	8 mm
 02-021S	02-021S SHORT NECK - STRAIGHT HOOK - LEFT	6,5 mm
 02-022S	02-022S SHORT NECK - STRAIGHT HOOK - RIGHT	6,5 mm
 02-023S	02-023S LONG NECK - STRAIGHT HOOK - LEFT	8 mm
 02-024S	02-024S LONG NECK - STRAIGHT HOOK - RIGHT	8 mm
 02-025S	02-025S SHORT NECK - LONG HOOK - LEFT	6,5 mm
 02-026S	02-026S SHORT NECK - LONG HOOK - RIGHT	6,5 mm
 02-027S	02-027S LONG NECK - LONG HOOK - LEFT	8 mm
 02-028S	02-028S LONG NECK - LONG HOOK - RIGHT	8 mm
 02-031S	02-031S PEDIATRIC - SHORT NECK - LONG HOOK - LEFT	6,5 mm
 02-032S	02-032S PEDIATRIC - SHORT NECK - LONG HOOK - RIGHT	6,5 mm



PEDIATRIC MODEL IS ONLY COMPATIBLE WITH 1.8 mm SCREWS

SCALE 1:1



OBA MANDIBULA

TITANIUM GRADE 2

	REF.	DESCRIPTION	NECK
 02-007S	02-007S	LONG NECK - DOUBLE HOOK	8 mm
 02-041S	02-041S	LONG NECK - DOUBLE HOOK - LEFT	8 mm
 02-042S	02-042S	LONG NECK - DOUBLE HOOK - RIGHT	8 mm
 02-043S	02-043S	LONG NECK - STRAIGHT HOOK- LEFT	8 mm
 02-044S	02-044S	LONG NECK - STRAIGHT HOOK - RIGHT	8 mm
 02-045S	02-045S	LONG NECK - LONG HOOK - LEFT	8 mm
 02-046S	02-046S	LONG NECK - LONG HOOK - RIGHT	8 mm



INSTRUMENTS

OBA



REF.	DESCRIPTION
99-901A	HANDPIECE FOR ALL INSERTS



99-909S	INSERT - INTERNAL PENTAGON
---------	----------------------------



99-910S	INSERT - EXTERNAL PENTAGON
---------	----------------------------



99-915S	INSERT - SLOTTED
---------	------------------



DRILLS CAN BE FOUND ON PAGE 33 OF THE CATALOGUE



SURGI-TEC ORTHODONTIC BONE ANCHOR

SCREWS TO USE WITH OBA

SELF-DRILLING



REF.	DIMENSIONS	HEAD TYPE
70-505S	Ø 2.0 mm - LENGTH 5 mm	ROUND/SLOT
70-507S	Ø 2.0 mm - LENGTH 7 mm	ROUND/SLOT



SELF-TAPPING



70-105S	Ø 2.3 mm - LENGTH 5 mm	ROUND/SLOT
70-107S	Ø 2.3 mm - LENGTH 7 mm	ROUND/SLOT



70-205S	Ø 2.3 mm - LENGTH 5 mm	LOW PENTAGON
70-207S	Ø 2.3 mm - LENGTH 7 mm	LOW PENTAGON



SELF-TAPPING "PEDIATRIC" (ONLY FOR USE WITH 02-03 IS AND 02-032S)




73-305S	Ø 1.8 mm - LENGTH 5 mm	ROUND/SLOT
73-307S	Ø 1.8 mm - LENGTH 7 mm	ROUND/SLOT



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5

OSTEOSYNTHESIS PLATES



SLA: SAND
BLASTED &
ACID ETCHED

0.6 mm PLATE THICKNESS

UNIQUE SLOT DESIGN
ALLOWS FINE ADJUSTMENT

Ø 2.4 mm

5.4 mm MOVABLE SLOTS

Ø 1.4 mm

SEMI-RIGID TITANIUM GRADE I
ROUGH SURFACE PREVENTS SLIPPING AND ALLOWS
FOR FASTER HEALING

NO SHARP EDGES

Ø 2.4 mm

EXTRA LOCKING
AFTER PERFECT PLACEMENT

SURGI-TEC 

THE ART OF MOVING BONE



SURGI-TEC OSTEOSYNTHESIS PLATES

OSP

TITANIUM GRADE I

	REF.	DESCRIPTION	TOTAL HEIGHT
	60-001S	SLOT PLATE - SLA - DIAMOND	17.60 mm
	60-003S	SLOT PLATE - SLA - RUBY	21.05 mm
	60-004S	SLOT PLATE - SLA - EMERALD	26.05 mm
	60-005S	SLOT PLATE - SLA - TOPAZ	24.00 mm
	60-006S	SLOT PLATE - SLA - AMETHYST	30.00 mm
	60-007S	SLOT PLATE - SLA - GENIOPLASTY - OPAL	14.60 mm
	60-008S	SLOT PLATE - SLA - GENIOPLASTY - ONYX	18.10 mm



INSTRUMENTS

OSP

TO USE WITH 1.8 mm, 2.0 mm, 2.3 mm & 2.5 mm SCREWS

	REF.	DESCRIPTION
	99-901A	HANDPIECE FOR ALL INSERTS
	99-909S	INSERT - INTERNAL PENTAGON
	99-910S	INSERT - EXTERNAL PENTAGON
	99-101A	SCREWDRIVER INTERNAL PENTAGON
	99-915S	INSERT - SLOTTED

TO USE WITH 1.2 mm & 1.4 mm MICRO SCREWS

	99-901A	HANDPIECE FOR ALL INSERTS
	99-906S	INSERT - MICRO INTERNAL PENTAGON
	99-906A	SCREWDRIVER - MICRO INTERNAL PENTAGON
	99-907A	SCREWDRIVER - MICRO EXTERNAL PENTAGON



DRILLS CAN BE FOUND ON PAGE 33 OF THE CATALOGUE

SCREWS TO USE WITH OSP

SELF-TAPPING



REF.	DIMENSIONS	HEAD TYPE
70-105S	Ø 2.3 mm - LENGTH 5 mm	ROUND/SLOT
70-107S	Ø 2.3 mm - LENGTH 7 mm	ROUND/SLOT



70-805S	Ø 2.3 mm - LENGTH 5 mm	PENTAGON
70-807S	Ø 2.3 mm - LENGTH 7 mm	PENTAGON
70-809S	Ø 2.3 mm - LENGTH 9 mm	PENTAGON
70-811S	Ø 2.3 mm - LENGTH 11 mm	PENTAGON
70-813S	Ø 2.3 mm - LENGTH 13 mm	PENTAGON
70-815S	Ø 2.3 mm - LENGTH 15 mm	PENTAGON
70-817S	Ø 2.3 mm - LENGTH 17 mm	PENTAGON
70-819S	Ø 2.3 mm - LENGTH 19 mm	PENTAGON

MICRO SCREW



71-205S	Ø 1.2 mm - LENGTH 5 mm	PENTAGON
71-207S	Ø 1.2 mm - LENGTH 7 mm	PENTAGON

RESCUE SCREW




70-905S	Ø 2.5 mm - LENGTH 5 mm	PENTAGON
70-907S	Ø 2.5 mm - LENGTH 7 mm	PENTAGON



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides.

If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5

OSTEOSYNTHESIS SCREWS





SURGI-TEC OSTEOSYNTHESIS SCREWS

OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

SELF-DRILLING



REF.	DIMENSIONS	HEAD TYPE
------	------------	-----------

70-505S Ø 2.0 mm - LENGTH 5 mm ROUND/SLOT



70-507S Ø 2.0 mm - LENGTH 7 mm ROUND/SLOT



70-707S Ø 2.0 mm - LENGTH 7 mm PENTAGON



70-709S Ø 2.0 mm - LENGTH 9 mm PENTAGON

SELF-TAPPING



70-605S Ø 2.0 mm - LENGTH 5 mm ROUND/SLOT



70-607S Ø 2.0 mm - LENGTH 7 mm ROUND/SLOT



70-609S Ø 2.0 mm - LENGTH 9 mm ROUND/SLOT



70-611S Ø 2.0 mm - LENGTH 11 mm ROUND/SLOT



70-105S Ø 2.3 mm - LENGTH 5 mm ROUND/SLOT



70-107S Ø 2.3 mm - LENGTH 7 mm ROUND/SLOT



70-205S Ø 2.3 mm - LENGTH 5 mm LOW PENTAGON



70-207S Ø 2.3 mm - LENGTH 7 mm LOW PENTAGON



70-217S Ø 2.3 mm - LENGTH 17 mm LOW PENTAGON



70-219S Ø 2.3 mm - LENGTH 19 mm LOW PENTAGON



70-221S Ø 2.3 mm - LENGTH 21 mm LOW PENTAGON



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5

OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

SELF-TAPPING

	REF.	DIMENSIONS	HEAD TYPE
 	70-805S	Ø 2.3 mm - LENGTH 5 mm	PENTAGON
 	70-807S	Ø 2.3 mm - LENGTH 7 mm	PENTAGON
 	70-809S	Ø 2.3 mm - LENGTH 9 mm	PENTAGON
 	70-811S	Ø 2.3 mm - LENGTH 11 mm	PENTAGON
 	70-813S	Ø 2.3 mm - LENGTH 13 mm	PENTAGON
 	70-815S	Ø 2.3 mm - LENGTH 15 mm	PENTAGON
 	70-817S	Ø 2.3 mm - LENGTH 17 mm	PENTAGON
 	70-819S	Ø 2.3 mm - LENGTH 19 mm	PENTAGON

SELF-TAPPING "RESCUE"

 	70-905S	Ø 2.5 mm - LENGTH 5 mm	PENTAGON
 	70-907S	Ø 2.5 mm - LENGTH 7 mm	PENTAGON

SELF-DRILLING "MICRO-SUSPENSION"

 	71-305S	Ø 1.4 mm - LENGTH 5 mm	PENTAGON
 	71-307S	Ø 1.4 mm - LENGTH 7 mm	PENTAGON



Microsuspensions screws are to be used with a micro internal pentagon screwdriver (99-906S)

SELF-TAPPING "PEDIATRIC"

 	73-305S	Ø 1.8 mm - LENGTH 5 mm	ROUND/SLOT
 	73-307S	Ø 1.8 mm - LENGTH 7 mm	ROUND/SLOT



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5



SURGI-TEC OSTEOSYNTHESIS SCREWS


OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

	REF.	DIMENSIONS	HEAD TYPE
SELF-TAPPING MICRO-SCREWS FOR “BONE GRAFTS”			
 	71-205S	Ø 1.2 mm - LENGTH 5 mm	PENTAGON
 	71-207S	Ø 1.2 mm - LENGTH 7 mm	PENTAGON
 	71-209S	Ø 1.2 mm - LENGTH 9 mm	PENTAGON
 	71-211S	Ø 1.2 mm - LENGTH 11 mm	PENTAGON
 	71-213S	Ø 1.2 mm - LENGTH 13 mm	PENTAGON
 	71-215S	Ø 1.2 mm - LENGTH 15 mm	PENTAGON
 	71-217S	Ø 1.2 mm - LENGTH 17 mm	PENTAGON



Microsuspensions screws are to be used with a micro internal pentagon screwdriver (99-906S)



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5



SURGI-TEC OSTEOSYNTHESIS DRILLS

DRILLS

CUSTOMIZED FOR SURGI-TEC PRODUCTS

2 cutting edges
stainless steel



REF.	DIAMETER	USABLE LENGTH	TOTAL LENGTH
46050/07/010	1.00 mm	7 mm	50 mm
46083/22/010	1.00 mm	22 mm	83 mm
46050/15/012	1.20 mm	15 mm	50 mm
46050/07/0165	1.65 mm	7 mm	50 mm
46098/20/0165	1.65 mm	20 mm	98 mm
46050/07/018	1.80 mm	7 mm	50 mm
46098/25/018	1.80 mm	25 mm	98 mm

SURGI-TEC'S COLOUR CODING

Ø 1.00 mm	RED	1 RING	FOR Ø 1.2 mm SCREWS
Ø 1.20 mm	RED	2 RINGS	FOR DRILLING A GLIDING HOLE Ø 1.2 mm SCREWS (LAG SCREW)
Ø 1.65 mm	BLUE	1 RING	FOR Ø 2.3 mm SCREWS IN THE MAXILLA
Ø 1.80 mm	BROWN	1 RING	FOR Ø 2.3 mm SCREWS IN THE MANDIBLE



Surgi-tec Recommends the following Ødrills for its different screws:
 Ø1mm drills: are commonly used with the micro screws (Ø1.2mm)
 Ø1.65mm drills: are recommended with Ø2.3mm screws in the Maxilla
 where Ø1.8mm drills are recommended with Ø2.3mm screws in the Mandible.
 We offer slight deviations on these recommended drills to suit all.

Note: we trust the surgeon's knowledge to choose a suitable drill size.

MOMMAERTS ORTHOFACIAL INSTRUMENTS

REF.	DESCRIPTION
------	-------------

MOMMAERTS LE FORT I RETRACTOR



99-201W

THE ANATOMICAL SHAPE OF THE RETRACTOR'S NECK ALLOWS COMFORTABLE ASSISTANCE OF THE HANDLING OF A MOBILISED MAXILLA IN A LE FORT I OSTEOTOMY. A SLIPPING RETRACTOR WILL NOT HURT THE LOWER LIP, BECAUSE THE TIP OF THE HOOK IS BLUNT.

MOMMAERTS - ABELOOS CHIN RETRACTOR



99-203W

CHIN SEGMENT HOLDER AND LIP RETRACTION WITHOUT UNDUE DEGLOVING. STABLE AND COMFORTABLE ASSISTANCE DURING A SLIDING GENIOPLASTY.

MOMMAERTS ORTHOFACIAL INSTRUMENTS

REF.

DESCRIPTION

MOMMAERTS ENDO BSSO RETRACTOR



99-202W

SUITED FOR THE LINGUAL AND BUCCAL TUNNEL OF ENDOSCOPICALLY ASSISTED, BUT ALSO OPEN SKY SAGITTAL SPLIT PROCEDURES. LIMITED DISSECTION AND EXCELLENT PROTECTION OF THE PERIOSTEUM FROM ROTATING INSTRUMENTS.

MOMMAERTS PERIOSTUM SCRAPER



99-204W

THE BAYONET SIDE IS DESIGNED FOR THE LINGUAL TUNNEL PREPARATION OF ENDOSCOPICALLY ASSISTED OR OPEN SKY SAGITTAL SPLIT PROCEDURES. THE OTHER SIDE IS MEANT TO PROTECT THE MARGINAL BRANCH OF THE FACIAL NERVE DURING SLIDING GENIOPLASTY PROCEDURES.

SURGI - TEC INSTRUMENTS

	REF.	DESCRIPTION
	03-750S	TPD NEO PATIENT KEY
	03-950S	TPD ALL-IN-ONE / CLASSIC PATIENT KEY
	03-751A	TPD NEO HINGED KEY
	03-951A	TPD ALL-IN-ONE / CLASSIC HINGED KEY

SURGI - TEC INSTRUMENTS

	REF.	DESCRIPTION
	06-950A	TMD PATIENT KEY
	99-901A	HANDPIECE FOR ALL INSERTS
	99-909S	INSERT - INTERNAL PENTAGON
	99-101A	SCREWDRIVER - INTERNAL PENTAGON
	99-910S	INSERT - EXTERNAL PENTAGON
	99-915S	INSERT - SLOTTED
	99-906S	INSERT - MICRO INTERNAL PENTAGON
	99-906A	SCREWDRIVER - MICRO INTERNAL PENTAGON
	99-907A	SCREWDRIVER - MICRO EXTERNAL PENTAGON

medartis® EXCLUSIVE DISTRIBUTOR

PRECISION IN FIXATION

As an innovation leader, Medartis promotes technologies and solutions for osteosynthesis in the field of cranio-maxillofacial surgery as well as upper and lower extremities.

We are the exclusive distributor of the Medartis MODUS products in Belgium.

Please scan the QR code to check their catalogue.





The simple yet comprehensive orthodontic mini implant system for an infinite number of anchorage applications. DB Orthodontics strives for absolute perfection from the design stage to the manufacturing of the Infinitas Mini Implant system.

Please scan the QR code to check their catalogue.





We are a distributor of META, the renowned company that produces sterile, single use bone harvesting instruments.

Please scan the QR code to check their catalogue.



META has a conspicuous and active role in the medical and dental fields, thanks to the extreme care given to every single phase of its manufacturing activity. We are the distributor of their single use bone harvesting instruments: the micross, the safescraper and the new smartscraper.

These instruments provides an easy method to obtain cortical bone for grafts in any type of defects. It is a fast solution, cost-effective and ready to use as the disposable device is proposed in a sterile blister pack.



[illegible]

41

A photograph of a large, light-colored industrial building with a dark sign that reads "SURGI-TEC" and a circular logo. The building is partially obscured by a dense green hedge in the foreground. A large teal triangle is on the left side of the image, and the sky is blue with white clouds.

SURGI-TEC

THE ART OF MOVING BONE



WWW.SURGI-TEC.COM

ALL RIGHTS RESERVED TO SURGI-TEC NV



MedTech
Flanders



Flemish
Government



CADskills



Flemish
Institute for
Entrepreneurship