

SURGI-TEC

THE ART OF MOVING BONE



WWW.SURGI-TEC.COM

SURGI-TEC

THE ART OF MOVING BONE



Surgi-Tec is an ISO 13485 certified medical device manufacturer based in Belgium. We specialise in CMF surgery implants.

The development of the TPD (transpalatal distractor) and the OBA (ortho bone anchor) have caused a remarkable shift in the field of CMF surgery.

These products are developed in close collaboration with Prof. Dr. M.Y. Mommaerts. This made the company renowned in Belgium and its neighbouring countries.

Quality is our top priority: we only work with certified suppliers for all our materials.

Our high-precision machines and the expertise of our production team assures that these quality standards are guaranteed in our assembly process.

The same quality requirements are continued in the relationship with our partners and customer service.


Personal contact is one of the core values of Surgi-Tec.

Our innovation and know-how have led to a very extensive OEM and subcontracting business.

Your idea is our common project!

SURGI-TEC NV BELGIUM
POORTAKKERSTRAAT 43-9051 SINT-DENIJS-WESTREM
PHONE: +32 (0)9 230 05 69
WWW.SURGI-TEC.COM - INFO@SURGI-TEC.COM

OSTEOSYNTHESIS SCREWS



SURGI-TEC 

THE ART OF MOVING BONE



OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

SELF-DRILLING



REF.	DIMENSIONS	HEAD TYPE
------	------------	-----------

70-505S	Ø 2.0 mm - LENGTH 5 mm	ROUND/SLOT
---------	------------------------	------------



70-507S	Ø 2.0 mm - LENGTH 7 mm	ROUND/SLOT
---------	------------------------	------------



70-707S	Ø 2.0 mm - LENGTH 7 mm	PENTAGON
---------	------------------------	----------



70-709S	Ø 2.0 mm - LENGTH 9 mm	PENTAGON
---------	------------------------	----------

SELF-TAPPING



70-605S	Ø 2.0 mm - LENGTH 5 mm	ROUND/SLOT
---------	------------------------	------------



70-607S	Ø 2.0 mm - LENGTH 7 mm	ROUND/SLOT
---------	------------------------	------------



70-609S	Ø 2.0 mm - LENGTH 9 mm	ROUND/SLOT
---------	------------------------	------------



70-611S	Ø 2.0 mm - LENGTH 11 mm	ROUND/SLOT
---------	-------------------------	------------



70-105S	Ø 2.3 mm - LENGTH 5 mm	ROUND/SLOT
---------	------------------------	------------



70-107S	Ø 2.3 mm - LENGTH 7 mm	ROUND/SLOT
---------	------------------------	------------



70-205S	Ø 2.3 mm - LENGTH 5 mm	LOW PENTAGON
---------	------------------------	--------------



70-207S	Ø 2.3 mm - LENGTH 7 mm	LOW PENTAGON
---------	------------------------	--------------



70-217S	Ø 2.3 mm - LENGTH 17 mm	LOW PENTAGON
---------	-------------------------	--------------



70-219S	Ø 2.3 mm - LENGTH 19 mm	LOW PENTAGON
---------	-------------------------	--------------



70-221S	Ø 2.3 mm - LENGTH 21 mm	LOW PENTAGON
---------	-------------------------	--------------



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5



OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

SELF-TAPPING



REF.	DIMENSIONS	HEAD TYPE
70-805S	Ø 2.3 mm - LENGTH 5 mm	PENTAGON
70-807S	Ø 2.3 mm - LENGTH 7 mm	PENTAGON
70-809S	Ø 2.3 mm - LENGTH 9 mm	PENTAGON
70-811S	Ø 2.3 mm - LENGTH 11 mm	PENTAGON
70-813S	Ø 2.3 mm - LENGTH 13 mm	PENTAGON
70-815S	Ø 2.3 mm - LENGTH 15 mm	PENTAGON
70-817S	Ø 2.3 mm - LENGTH 17 mm	PENTAGON
70-819S	Ø 2.3 mm - LENGTH 19 mm	PENTAGON



SELF-TAPPING "RESCUE"



70-905S	Ø 2.5 mm - LENGTH 5 mm	PENTAGON
70-907S	Ø 2.5 mm - LENGTH 7 mm	PENTAGON



SELF-DRILLING "MICRO-SUSPENSION"



71-305S	Ø 1.4 mm - LENGTH 5 mm	PENTAGON
71-307S	Ø 1.4 mm - LENGTH 7 mm	PENTAGON



Microsuspensions screws are to be used with a micro internal pentagon screwdriver (99-906S)

SELF-TAPPING "PEDIATRIC"



73-305S	Ø 1.8 mm - LENGTH 5 mm	ROUND/SLOT
73-307S	Ø 1.8 mm - LENGTH 7 mm	ROUND/SLOT



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our sets and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5



SURGI-TEC

THE ART OF MOVING BONE



WWW.SURGI-TEC.COM

ALL RIGHTS RESERVED TO SURGI-TEC NV

