

SURGI-TEC



THE ART OF MOVING BONE



WWW.SURGI-TEC.COM

SURGI-TEC

THE ART OF MOVING BONE



Surgi-Tec is a Belgian company with a reach into the international medical industry. More than twenty years ago, the company was established with a clear mission: to produce the best medical devices in titanium for maxillofacial surgery, and improve the life quality of patients worldwide.

Since then, Surgi-Tec mastered the art of moving bone and developed life-changing innovations such as the transpalatal distractor and orthodontic bone anchors.

Our surgical instruments and devices are top-notch, and the same quality standards are continued in the relationship with our partners and customers. Personal contact is one of Surgi-Tec's core values.

Surgi-Tec works only with certified suppliers for all materials. Our high-precision machinery and the expertise of our production team ensure that the highest quality standards are guaranteed in our assembly process. We also supply our implantable medical devices in sterile packaging. Please contact the Surgi-Tec sales team if required.

Our mentality, together with our innovation and know-how have also led to a very extensive OEM and subcontracting business: Your idea is our joint project!

SURGI-TEC NV BELGIUM
POORTAKKERSTRAAT 43-9051 SINT-DENIJS-WESTREM
PHONE: +32 (0)9 230 05 69
WWW.SURGI-TEC.COM - INFO@SURGI-TEC.COM

TABLE OF CONTENTS

TRANSPALATAL DISTRACTORS

| | |
|----------------------------|----|
| TPD NEO | 7 |
| TPD ALL-IN-ONE | 9 |
| TPD CLASSIC | 11 |
| INSTRUMENTS AND SCREWS TPD | 12 |

ORTHODONTIC BONE ANCHORS

| | |
|----------------------------|----|
| MAXILLA | 16 |
| MANDIBLE AND OBA SLA | 17 |
| INSTRUMENTS AND SCREWS OBA | 18 |

OSTEOSYNTHESIS PLATES

| | |
|----------------------------|----|
| OSP | 21 |
| INSTRUMENTS AND SCREWS OSP | 22 |

OSTEOSYNTHESIS SCREWS 26


INSTRUMENTS

| | |
|-------------------------|----|
| DRILLS | 29 |
| ORTHOFACIAL INSTRUMENTS | 30 |

DISTRIBUTION

| | |
|------------------------------|----|
| TTS - TITANIUM TRAUMA SPLINT | 32 |
| INFINITAS | 33 |
| META | 34 |

TRANS PALATAL DISTRACTOR NEO



EXTRA SHARP
SPIKES

EASY PLACEMENT

COUNTER NUT

CLOSED BODY DESIGN

COLOR MARKING
1 STEP IS 0.33 mm DISTRACTION



3 DIMENSIONAL MOVEMENT

INDICATIONS FOR
PLACEMENT

ROUNDED EDGES

SLIM AND CLOSED DESIGN INCREASES PATIENT COMFORT

SURGI-TEC 





THE ART OF MOVING BONE



SURGI-TEC TRANSPALATAL DISTRACTOR

TPD NEO

TITANIUM GRADE 2 & 5

| | REF. | DESCRIPTION | INITIAL WIDTH | DISTRACTION | MAX. WIDTH |
|---|---------|-------------------------|---------------|-------------|------------|
|  | 03-701A | TPD NEO - MODULE 1 | 17 mm | 7 mm | 24 mm |
|  | 03-702A | TPD NEO - MODULE 2 | 20 mm | 12 mm | 32 mm |
|  | 03-725A | TPD NEO - MODULE 2.5 | 24 mm | 19 mm | 43 mm |
|  | 03-703A | TPD NEO - MODULE 3 | 27 mm | 25 mm | 52 mm |


SCALE 1:1



Surgi-Tec TPD NEO is also available in sterile packaging.
If you prefer to receive the device in sterile packaging, please add "- STER" to the reference.
e.g.: 03-725A in sterile packaging becomes 03-725A - STER

TRANS PALATAL DISTRACTOR

ALL-IN-ONE



EXTRA SHARP
SPIKES

LOCKING SCREW



INDICATIONS
FOR PLACEMENT

COLOR MARKING
1 STEP IS 0.33 mm DISTRACTION

ROUNDED EDGES

3 DIMENSIONAL MOVEMENT



TELESCOPIC DESIGN FOR
MAXIMUM DISTRACTION

SURGI-TEC

THE ART OF MOVING BONE




SURGI-TEC TRANSPALATAL DISTRACTOR

TPD ALL-IN-ONE

TITANIUM GRADE 2 & 5

| | REF. | DESCRIPTION | INITIAL WIDTH | DISTRACTION | MAX. WIDTH |
|---|---------|---|-------------------|-------------|------------|
|  | 03-901A | TPD All-In-One - MODULE 1 | 15 mm | 7 mm | 22 mm |
|  | 03-902A | TPD All-In-One - MODULE 2 | 18 mm | 12 mm | 30 mm |
|  | 03-925A | TPD All-In-One - MODULE 2.5 | 21 mm | 18 mm | 39 mm |
|  | 03-903A | TPD All-In-One - MODULE 3 | 24 mm | 23 mm | 47 mm |
|  | 03-904A | TPD All-In-One - MODULE 4 | 29 mm | 29 mm | 58 mm |
|  | 99-100S | TPD - All In One / Classic LOCKING SCREW | INTERNAL PENTAGON | | |

TRANS PALATAL DISTRACTOR CLASSIC



EXTRA SHARP
SPIKES

LOCKING SCREW



SMOOTH SURFACE

COLOR MARKING
1 STEP IS 0.33 mm DISTRACTION

ROUNDED EDGES

HOLES FOR LIGATURE

3 PIECE DESIGN ALLOWS FOR INTERCHANGEABLE MODULES

SURGI-TEC 







THE ART OF MOVING BONE



SURGI-TEC TRANSPALATAL DISTRACTOR

TPD CLASSIC

TITANIUM GRADE 2 & 5



| | REF. | DESCRIPTION | INITIAL WIDTH | DISTRACTION | MAX. WIDTH |
|---|---------|--|---------------|-------------------|------------|
|  | 03-800A | TPD CLASSIC - LEFT AND RIGHT ABUTMENT PLATES | | 2.5 mm holes | |
|  | 03-801A | TPD CLASSIC - MODULE 1 | 15 mm | 7 mm | 22 mm |
|  | 03-802A | TPD CLASSIC - MODULE 2 | 18 mm | 12 mm | 30 mm |
|  | 03-825A | TPD CLASSIC - MODULE 2.5 | 21 mm | 18 mm | 39 mm |
|  | 03-803A | TPD CLASSIC - MODULE 3 | 24 mm | 23 mm | 47 mm |
|  | 99-100S | TPD - All In One / Classic LOCKING SCREW | | INTERNAL PENTAGON | |



SURGI-TEC TRANS PALATAL DISTRACTOR




SURGICAL INSTRUMENTS

TPD NEO

| | REF. | DESCRIPTION |
|---|---------|----------------------------|
|  | 03-751A | TPD NEO - HINGED KEY |
|  | 99-901A | HANDPIECE FOR ALL INSERTS |
|  | 99-909S | INSERT - INTERNAL PENTAGON |
|  | 99-910S | INSERT - EXTERNAL PENTAGON |

INSTRUMENTS

TPD NEO

| | REF. | DESCRIPTION |
|---|---------|---|
|  | 03-750S | TPD NEO - PATIENT KEY |
|  | 03-750A | TPD NEO - HINGED PATIENT KEY |
|  | 99-701W | TPD NEO PLACING FORCEPS Optional accessory Surgical Sarl - distributed by SURGI-TEC |



SURGI-TEC TRANS PALATAL DISTRACTOR

SURGICAL INSTRUMENTS

TPD ALL-IN-ONE AND TPD CLASSIC

| | REF. | DESCRIPTION |
|---|---------|--|
|  | 03-951A | TPD ALL-IN-ONE / CLASSIC - HINGED KEY |
|  | 99-901A | HANDPIECE FOR ALL INSERTS |
|  | 99-101A | SCREWDRIVER - INTERNAL PENTAGON |
|  | 99-909S | INSERT - INTERNAL PENTAGON |
|  | 99-910S | INSERT - EXTERNAL PENTAGON |

INSTRUMENTS

TPD ALL-IN-ONE AND TPD CLASSIC

| | REF. | DESCRIPTION |
|---|---------|---|
|  | 03-950S | TPD ALL-IN-ONE / CLASSIC - PATIENT KEY |



SURGI-TEC TRANS PALATAL DISTRACTOR

INSTRUMENTS

DUMMIES (SCALE 1:1)






| | REF. | DESCRIPTION |
|--|---------|------------------------------------|
|  | 03-790S | TPD NEO - DUMMIES |
|  | 03-990S | TPD ALL-IN-ONE / CLASSIC - DUMMIES |

SCREWS

All TPD's are compatible with 2.5 mm rescue screws


TPD NEO and TPD All-In-One: 2 screws needed

TPD CLASSIC: 4 screws needed

| | REF. | DESCRIPTION |
|---|---------|--|
|   | 70-707S | Ø 2.0 Selfdrilling Screw - L 7 mm - Pentagon |
|   | 70-709S | Ø 2.0 Selfdrilling Screw - L 9 mm - Pentagon |
|   | 70-807S | Ø 2.3 Selftapping Screw - L 7 mm - Pentagon |

ORTHODONTIC BONE ANCHOR

EASILY ADAPTABLE TO THE BONE SURFACE



SURGI-TEC 









THE ART OF MOVING BONE



SURGI-TEC ORTHODONTIC BONE ANCHOR

OBA MAXILLA

TITANIUM GRADE 2

| | REF. | DESCRIPTION | NECK |
|--|---------|---|--------|
|  02-005S | 02-005S | OBA MAXILLA - LONG NECK - DOUBLE HOOK | 8 mm |
|  02-021S | 02-021S | OBA MAXILLA - SHORT NECK - STRAIGHT HOOK - LEFT | 6,5 mm |
|  02-022S | 02-022S | OBA MAXILLA - SHORT NECK - STRAIGHT HOOK - RIGHT | 6,5 mm |
| 02-021S 02-022S | | | |
|  02-023S | 02-023S | OBA MAXILLA - LONG NECK - STRAIGHT HOOK - LEFT | 8 mm |
|  02-024S | 02-024S | OBA MAXILLA - LONG NECK - STRAIGHT HOOK - RIGHT | 8 mm |
| 02-023S 02-024S | | | |
|  02-025S | 02-025S | OBA MAXILLA - SHORT NECK - LONG HOOK - LEFT | 6,5 mm |
|  02-026S | 02-026S | OBA MAXILLA - SHORT NECK - LONG HOOK - RIGHT | 6,5 mm |
| 02-025S 02-026S | | | |
|  02-027S | 02-027S | OBA MAXILLA - LONG NECK - LONG HOOK - LEFT | 8 mm |
|  02-028S | 02-028S | OBA MAXILLA - LONG NECK - LONG HOOK - RIGHT | 8 mm |
| 02-027S 02-028S | | | |
|  02-031S | 02-031S | OBA MAXILLA - PEDIATRIC - SHORT NECK - LONG HOOK - LEFT | 6,5 mm |
|  02-032S | 02-032S | OBA MAXILLA - PEDIATRIC - SHORT NECK - LONG HOOK - RIGHT | 6,5 mm |
| 02-031S 02-032S | | | |



PEDIATRIC MODEL IS ONLY COMPATIBLE WITH 1.8 mm SCREWS












SURGI-TEC ORTHODONTIC BONE ANCHOR





OBA MANDIBLE

TITANIUM GRADE 2

|  02-007S | REF. | DESCRIPTION | NECK |
|--|---------|--|------|
| | 02-007S | OBA MANDIBLE - LONG NECK - DOUBLE HOOK | 8 mm |
|   02-041S 02-042S | 02-041S | OBA MANDIBLE - LONG NECK - DOUBLE HOOK - LEFT | 8 mm |
| | 02-042S | OBA MANDIBLE - LONG NECK - DOUBLE HOOK - RIGHT | 8 mm |
|   02-043S 02-044S | 02-043S | OBA MANDIBLE - LONG NECK - STRAIGHT HOOK - LEFT | 8 mm |
| | 02-044S | OBA MANDIBLE - LONG NECK - STRAIGHT HOOK - RIGHT | 8 mm |
|   02-045S 02-046S | 02-045S | OBA MANDIBLE - LONG NECK - LONG HOOK - LEFT | 8 mm |
| | 02-046S | OBA MANDIBLE - LONG NECK - LONG HOOK - RIGHT | 8 mm |

OBA SLA - maxilla and mandible

TITANIUM GRADE 2 - SANDBLASTED AND ACID ETCHED

|   02-003S 02-006S | REF. | DESCRIPTION | NECK |
|--|---------|------------------------------------|---------|
| | 02-003S | OBA SLA - SHORT NECK - DOUBLE HOOK | 4,5 mm |
|   02-009S 02-050S | 02-006S | OBA SLA - LONG NECK - DOUBLE HOOK | 6,5 mm |
| | 02-009S | OBA SLA - ZYGOMA - DOUBLE HOOK | 12,5 mm |
| | 02-050S | OBA SLA - ZYGOMA - LONG HOOK | 14,5 mm |

SCALE 1:1








Surgi-Tec OBA's are also available in sterile packaging.
If you prefer to receive the device in sterile packaging, please add "- STER" to the reference.
e.g.: 02-043S in sterile packaging becomes 02-043S - STER



SURGI-TEC ORTHODONTIC BONE ANCHOR

SURGICAL INSTRUMENTS

OBA

| | REF. | DESCRIPTION |
|--|---------|----------------------------|
|  | 99-901A | HANDPIECE FOR ALL INSERTS |
|  | 99-909S | INSERT - INTERNAL PENTAGON |
|  | 99-915S | INSERT - SLOTTED |
|  | 99-910S | INSERT - EXTERNAL PENTAGON |
|  This screwdriver can only be used with pentagon screws | | |



DI-210702263 ORTHODONTIC BONE ANCHOR TRAY
NTOC - distributed by SURGI-TEC



















DRILLS CAN BE FOUND ON PAGE 29 OF THIS CATALOGUE




SURGI-TEC ORTHODONTIC BONE ANCHOR

SCREWS TO USE WITH OBA

| | | REF. | DESCRIPTION |
|---|---|---------|---|
| SELF-DRILLING | | | |
|  |  | 70-505S | Ø 2.0 Selfdrilling Screw - L 5 mm - Round |
|  |  | 70-507S | Ø 2.0 Selfdrilling Screw - L 7 mm - Round |
| SELF-TAPPING | | | |
|  |  | 70-105S | Ø 2.3 Selftapping Screw - L 5 mm - Round |
|  |  | 70-107S | Ø 2.3 Selftapping Screw - L 7 mm - Round |
|  |  | 70-205S | Ø 2.3 Selftapping Screw - L 5 mm - low Pentagon |
|  |  | 70-207S | Ø 2.3 Selftapping Screw - L 7 mm - low Pentagon |
| SELF-TAPPING "PEDIATRIC" (ONLY FOR USE WITH 02-03 IS AND 02-032S) | | | |
|  |  | 73-305S | Ø 1.8 Selftapping screw - L 5 mm - Round |
|  |  | 73-307S | Ø 1.8 Selftapping screw - L 7 mm - Round |



OSTEOSYNTHESIS PLATES



SURGI-TEC 








THE ART OF MOVING BONE



SURGI-TEC OSTEOSYNTHESIS PLATES

OSP

TITANIUM GRADE I

| | REF. | DESCRIPTION | TOTAL HEIGHT |
|---|---------|---------------------------------------|--------------|
|  | 60-001S | SLOT PLATE - SLA - DIAMOND | 17.60 mm |
|  | 60-003S | SLOT PLATE - SLA - RUBY | 21.05 mm |
|  | 60-004S | SLOT PLATE - SLA - EMERALD | 26.05 mm |
|  | 60-005S | SLOT PLATE - SLA - TOPAZ | 24.00 mm |
|  | 60-006S | SLOT PLATE - SLA - AMETHYST | 30.00 mm |
|  | 60-007S | SLOT PLATE - SLA - GENIOPLASTY - OPAL | 14.60 mm |
|  | 60-008S | SLOT PLATE - SLA - GENIOPLASTY - ONYX | 18.10 mm |

SCALE 1:1






INSTRUMENTS

OSP

TO USE WITH 1.8 mm, 2.0 mm, 2.3 mm & 2.5 mm SCREWS

| | REF. | DESCRIPTION |
|---|---------|-------------------------------|
|  | 99-901A | HANDPIECE FOR ALL INSERTS |
|  | 99-909S | INSERT - INTERNAL PENTAGON |
|  | 99-910S | INSERT - EXTERNAL PENTAGON |
|  | 99-101A | SCREWDRIVER INTERNAL PENTAGON |
|  | 99-915S | INSERT - SLOTTED |

TO USE WITH 1.2 mm & 1.4 mm MICRO SCREWS

| | | |
|---|---------|---------------------------------------|
|  | 99-906S | INSERT - MICRO INTERNAL PENTAGON |
|  | 99-906A | SCREWDRIVER - MICRO INTERNAL PENTAGON |
|  | 99-907A | SCREWDRIVER - MICRO EXTERNAL PENTAGON |



DRILLS CAN BE FOUND ON PAGE 29 OF THIS CATALOGUE

SCREWS TO USE WITH OSP

| REF. | DESCRIPTION |
|---|--|
| SELF-TAPPING | |
|   | 70-205S Ø 2.3 Selftapping Screw - L 5 mm - low Pentagon |
|   | 70-207S Ø 2.3 Selftapping Screw - L 7 mm - low Pentagon |
|   | 70-209S Ø 2.3 Selftapping Screw - L 9 mm - low Pentagon |
|   | 70-211S Ø 2.3 Selftapping Screw - L 11 mm - low Pentagon |
|   | 70-213S Ø 2.3 Selftapping Screw - L 13 mm - low Pentagon |
|   | 70-215S Ø 2.3 Selftapping Screw - L 15 mm - low Pentagon |
|   | 70-217S Ø 2.3 Selftapping Screw - L 17 mm - low Pentagon |
|   | 70-219S Ø 2.3 Selftapping Screw - L 19 mm - low Pentagon |
|   | 70-221S Ø 2.3 Selftapping Screw - L 21 mm - low Pentagon |
|   | 70-805S Ø 2.3 Selftapping Screw - L 5 mm - Pentagon |
|   | 70-807S Ø 2.3 Selftapping Screw - L 7 mm - Pentagon |
|   | 70-809S Ø 2.3 Selftapping Screw - L 9 mm - Pentagon |
|   | 70-811S Ø 2.3 Selftapping Screw - L 11 mm - Pentagon |
|   | 70-813S Ø 2.3 Selftapping Screw - L 13 mm - Pentagon |
|   | 70-815S Ø 2.3 Selftapping Screw - L 15 mm - Pentagon |
|   | 70-817S Ø 2.3 Selftapping Screw - L 17 mm - Pentagon |
|   | 70-819S Ø 2.3 Selftapping Screw - L 19 mm - Pentagon |



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our trays and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides.
If you prefer to receive our screws in 5-packed format, please add .5 to the reference.
e.g.: 70-205S in 5-pack becomes 70-205S.5




SURGI-TEC OSTEOSYNTHESIS PLATES

SCREWS TO USE WITH OSP

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

SELF-TAPPING



| | |
|---------|---|
| 75-203S | Ø 1.5 Selftapping Micro Screw - L 3 mm - Round |
| 75-205S | Ø 1.5 Selftapping Micro Screw - L 5 mm - Round |
| 75-207S | Ø 1.5 Selftapping Micro Screw - L 7 mm - Round |
| 75-209S | Ø 1.5 Selftapping Micro Screw - L 9 mm - Round |
| 75-211S | Ø 1.5 Selftapping Micro Screw - L 11 mm - Round |



| | |
|---------|--|
| 70-105S | Ø 2.3 Selftapping Screw - L 5 mm - Round |
| 70-107S | Ø 2.3 Selftapping Screw - L 7 mm - Round |

RESCUE SCREW



| | |
|---------|--|
| 70-905S | Ø 2.5 Selftapping rescue Screw - L 5 mm - Pentagon |
| 70-907S | Ø 2.5 Selftapping rescue Screw - L 7 mm - Pentagon |

MICRO SCREW




| | |
|---------|---|
| 71-205S | Ø 1.2 Selftapping Micro Screw - L 5 mm - Pentagon |
| 71-207S | Ø 1.2 Selftapping Micro Screw - L 7 mm - Pentagon |



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our trays and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides.
If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5

OSTEOSYNTHESIS SCREWS



All Surgi-Tec screws are available in single packed and 5-pack format. These containers are developed to fit in our trays and are colour marked. All critical information is clearly visible on top of the containers. The batch numbers are visible from the sides. If you prefer to receive our screws in 5-packed format, please add .5 to the reference. e.g.: 70-205S in 5-pack becomes 70-205S.5

SURGI-TEC

THE ART OF MOVING BONE



OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

SELF-DRILLING



| REF. | DESCRIPTION |
|---------|---|
| 70-505S | Ø 2.0 Selfdrilling Screw - L 5 mm - Round |
| 70-507S | Ø 2.0 Selfdrilling Screw - L 7 mm - Round |
| 70-707S | Ø 2.0 Selfdrilling Screw - L 7 mm - Pentagon |
| 70-709S | Ø 2.0 Selfdrilling Screw - L 9 mm - Pentagon |
| 71-305S | Ø 1.4 Selfdrilling Micro-suspension Screw - L 5 mm - Pentagon |
| 71-307S | Ø 1.4 Selfdrilling Micro-suspension Screw - L 7 mm - Pentagon |



Microsuspensions screws are to be used with a micro internal pentagon screwdriver (99-906A)







































SELF-TAPPING



| | |
|---------|---|
| 75-203S | Ø 1.5 Selftapping Micro Screw - L 3 mm - Round |
| 75-205S | Ø 1.5 Selftapping Micro Screw - L 5 mm - Round |
| 75-207S | Ø 1.5 Selftapping Micro Screw - L 7 mm - Round |
| 75-209S | Ø 1.5 Selftapping Micro Screw - L 9 mm - Round |
| 75-211S | Ø 1.5 Selftapping Micro Screw - L 11 mm - Round |
| 73-305S | Ø 1.8 Selftapping screw - L 5 mm - Round |
| 73-307S | Ø 1.8 Selftapping screw - L 7 mm - Round |
| 70-605S | Ø 2.0 Selftapping Screw - L 5 mm - Round |
| 70-607S | Ø 2.0 Selftapping Screw - L 7 mm - Round |
| 70-609S | Ø 2.0 Selftapping Screw - L 9 mm - Round |
| 70-611S | Ø 2.0 Selftapping Screw - L 11 mm - Round |
| 70-105S | Ø 2.3 Selftapping Screw - L 5 mm - Round |
| 70-107S | Ø 2.3 Selftapping Screw - L 7 mm - Round |

OSTEOSYNTHESIS SCREWS

TITANIUM GRADE 5

| | | REF. | DESCRIPTION |
|---|---|---------|--|
| SELF-TAPPING | | | |
|  |  | 70-205S | Ø 2.3 Selftapping Screw - L 5 mm - low Pentagon |
|  |  | 70-207S | Ø 2.3 Selftapping Screw - L 7 mm - low Pentagon |
|  |  | 70-209S | Ø 2.3 Selftapping Screw - L 9 mm - low Pentagon |
|  |  | 70-211S | Ø 2.3 Selftapping Screw - L 11 mm - low Pentagon |
|  |  | 70-213S | Ø 2.3 Selftapping Screw - L 13 mm - low Pentagon |
|  |  | 70-215S | Ø 2.3 Selftapping Screw - L 15 mm - low Pentagon |
|  |  | 70-217S | Ø 2.3 Selftapping Screw - L 17 mm - low Pentagon |
|  |  | 70-219S | Ø 2.3 Selftapping Screw - L 19 mm - low Pentagon |
|  |  | 70-221S | Ø 2.3 Selftapping Screw - L 21 mm - low Pentagon |
|  |  | 70-805S | Ø 2.3 Selftapping Screw - L 5 mm - Pentagon |
|  |  | 70-807S | Ø 2.3 Selftapping Screw - L 7 mm - Pentagon |
|  |  | 70-809S | Ø 2.3 Selftapping Screw - L 9 mm - Pentagon |
|  |  | 70-811S | Ø 2.3 Selftapping Screw - L 11 mm - Pentagon |
|  |  | 70-813S | Ø 2.3 Selftapping Screw - L 13 mm - Pentagon |
|  |  | 70-815S | Ø 2.3 Selftapping Screw - L 15 mm - Pentagon |
|  |  | 70-817S | Ø 2.3 Selftapping Screw - L 17 mm - Pentagon |
|  |  | 70-819S | Ø 2.3 Selftapping Screw - L 19 mm - Pentagon |
| SELF-TAPPING "RESCUE" | | | |
|  |  | 70-905S | Ø 2.5 Selftapping rescue Screw - L 5 mm - Pentagon |
|  |  | 70-907S | Ø 2.5 Selftapping rescue Screw - L 7 mm - Pentagon |

SCALE 2:1



MICRO-SCREWS

TITANIUM GRADE 5

REF.

DESCRIPTION

SELF-TAPPING MICRO-SCREWS FOR "BONE GRAFTS"



71-205S Ø 1.2 Selftapping Micro Screw - L 5 mm - Pentagon



71-207S Ø 1.2 Selftapping Micro Screw - L 7 mm - Pentagon



71-209S Ø 1.2 Selftapping Micro Screw - L 9 mm - Pentagon



71-211S Ø 1.2 Selftapping Micro Screw - L 11 mm - Pentagon



71-213S Ø 1.2 Selftapping Micro Screw - L 13 mm - Pentagon




71-215S Ø 1.2 Selftapping Micro Screw - L 15 mm - Pentagon



71-217S Ø 1.2 Selftapping Micro Screw - L 17 mm - Pentagon



Microsuspensions screws are to be used with a micro internal pentagon screwdriver (99-906S)






SURGI-TEC OSTEOSYNTHESIS DRILLS

DRILLS

Risa GmbH - distributed by SURGI-TEC

CUSTOMIZED FOR SURGI-TEC PRODUCTS

2 cutting edges
stainless steel



2/1

1/1

| REF. | DIAMETER | USABLE LENGTH | TOTAL LENGTH |
|---------------|----------|---------------|--------------|
| 46050/07/010 | 1.00 mm | 7 mm | 50 mm |
| 46083/22/010 | 1.00 mm | 22 mm | 83 mm |
| 46050/15/012 | 1.20 mm | 15 mm | 50 mm |
| 46050/07/0165 | 1.65 mm | 7 mm | 50 mm |
| 46098/20/0165 | 1.65 mm | 20 mm | 98 mm |
| 46050/07/018 | 1.80 mm | 7 mm | 50 mm |
| 46098/25/018 | 1.80 mm | 25 mm | 98 mm |

Recommendations for use

Ø 1.00 mm
Ø 1.20 mm
Ø 1.65 mm
Ø 1.80 mm

FOR Ø 1.2 mm SCREWS
FOR Ø 1.2 mm SCREWS - GLIDING HOLE (LAG SCREW)
FOR Ø 2.3 mm SCREWS IN THE MAXILLA
FOR Ø 2.3 mm SCREWS IN THE MANDIBLE



Note that the operating surgeon is always responsible for choosing the appropriate drill size for each specific case

We offer slight deviations on these recommended drills to suit all.

ORTHOFACIAL INSTRUMENTS

MD Medical - distributed by SURGI-TEC

REF.

DESCRIPTION

MOMMAERTS LE FORT I RETRACTOR



99-201W

THE ANATOMICAL SHAPE OF THE RETRACTOR'S NECK ALLOWS COMFORTABLE ASSISTANCE OF THE HANDLING OF A MOBILISED MAXILLA IN A LE FORT I OSTEOTOMY. A SLIPPING RETRACTOR WILL NOT HURT THE LOWER LIP, BECAUSE THE TIP OF THE HOOK IS BLUNT.

MOMMAERTS - ABELOOS CHIN RETRACTOR



99-203W

CHIN SEGMENT HOLDER AND LIP RETRACTION WITHOUT UNDUE DEGLOVING. STABLE AND COMFORTABLE ASSISTANCE DURING A SLIDING GENIOPLASTY.

ORTHOFACIAL INSTRUMENTS

MD Medical - distributed by SURGI-TEC

REF.

DESCRIPTION

MOMMAERTS ENDO BSSO RETRACTOR



99-202W

SUITED FOR THE LINGUAL AND BUCCAL TUNNEL OF ENDOSCOPICALLY ASSISTED, BUT ALSO OPEN SKY SAGITTAL SPLIT PROCEDURES. LIMITED DISSECTION AND EXCELLENT PROTECTION OF THE PERIOSTEUM FROM ROTATING INSTRUMENTS.

MOMMAERTS PERIOSTUM SCRAPER



99-204W

THE BAYONET SIDE IS DESIGNED FOR THE LINGUAL TUNNEL PREPARATION OF ENDOSCOPICALLY ASSISTED OR OPEN SKY SAGITTAL SPLIT PROCEDURES. THE OTHER SIDE IS MEANT TO PROTECT THE MARGINAL BRANCH OF THE FACIAL NERVE DURING SLIDING GENIOPLASTY PROCEDURES.



medartis®

PRECISION IN FIXATION

TTS – Titanium Trauma Splint



An ideal method for splinting of traumatized teeth



As an innovation leader, Medartis promotes technologies and solutions for osteosynthesis in the field of cranio-maxillofacial surgery as well as upper and lower extremities.

We are the exclusive distributor of the Medartis TTS - Titanium Trauma Splints in Belgium.

Please scan the QR code to check the surgical technique.



| Art. No. | Description | Color | Pieces/ Pkg |
|----------|--------------|------------|-------------|
| TTS-1002 | 52 × 0.2 mm | matt, gold | 1 |
| TTS-1004 | 100 × 0.2 mm | matt, gold | 1 |

Material: Titanium (ASTM F67)
Thickness: 0.2 mm

EXCLUSIVE DISTRIBUTOR




SURGI-TEC DISTRIBUTION IN BELGIUM



The simple yet comprehensive orthodontic mini implant system for an infinite number of anchorage applications.

DB Orthodontics strives for absolute perfection from the design stage to the manufacturing of the Infinitas Mini Implant system.

Please scan the QR code to check their catalogue.





META

ra

We are a distributor of META, the renowned company that produces sterile, single use bone harvesting instruments.

Please scan the QR code to check their catalogue.



META has a conspicuous and active role in the medical and dental fields, thanks to the extreme care given to every single phase of its manufacturing activity. We are the distributor of their single use bone harvesting instruments: micross and safescraper.

These instruments provides an easy method to obtain cortical bone for grafts in any type of defects.

It is a fast solution, cost-effective and ready to use as the disposable device is proposed in a sterile blister pack.



A photograph of a large, light-colored industrial building with a dark sign that reads "SURGI-TEC" and a circular logo. In the foreground, there is a dense green hedge. The sky is blue with white clouds. A large teal triangle is on the left side of the image.

SURGI-TEC

THE ART OF MOVING BONE

WWW.SURGI-TEC.COM

ALL RIGHTS RESERVED TO SURGI-TEC NV